


PLANMECA PM2000 UPHOLSTERY KITS

We do not recover your old kit – we make you a new kit!
No Core Return Required

FINESSE
<ul style="list-style-type: none"> ➤ For Comfort ➤ Has Slight Detailing

<ul style="list-style-type: none"> ❑ All Vinyl Outer Cover ❑ All New High Quality Foam ❑ Stitching on Face of Back and Seat ❑ Clear Toe Cover

PLUSH
<ul style="list-style-type: none"> ➤ For ULTRA COMFORT During Long Procedures ➤ More Luxurious

<ul style="list-style-type: none"> ❑ All Vinyl Outer Cover PLUS Upgrade Naugasoft Vinyl ❑ All New High Quality Foam PLUS Upgrade Memory Foam ❑ Stitching on Face of Back and Seat ❑ Clear Toe Cover

SEAMLESS
<ul style="list-style-type: none"> ➤ For Comfort, and ➤ Easy Cleaning

<ul style="list-style-type: none"> ❑ All Vinyl Outer Cover ❑ All New High Quality Foam ❑ No Stitching on Face of Back and Seat ❑ Clear Toe Cover



“OUR MOST POPULAR KIT”

The **PLUSH** style was developed at the request of many doctors seeking more patient comfort during long procedures.

Aerostatic achieves this comfort level by using a Super Soft **Memory Foam**. The heat of ones body causes the memory foam to mold to the body, giving much more uniform support, thus more comfort. In addition, Aerostatic uses the upgraded softer vinyls which allow the memory foam to function as designed.

TO ORDER YOUR KIT PACKAGE, PLEASE PROVIDE:

STYLE CHOSEN:

VINYL CHOSEN:

AEROSTATIC TECHNOLOGY CANADA

Phone: 1-888-891-1052 Fax: (519) 833-1053

E-mail: aerostatic@rogers.com

VINYL QUALITY!

It should not be assumed that the quality of all vinyls used in the Dental industry are the same. One cannot make an intelligent choice merely by thickness and feel of the product alone. Some vinyls have characteristics that are more beneficial in a Dental office environment. The following chart outlines the performance tests on the vinyls that Aerostatic recommends.

QUALITY AND PRICE COMPARISON FOR VINYL

PRODUCT	PRICE	TREATMENT TO RESIST						Results Abrasion Test
		Bacteria	Microbial	** Fungal	Sulfide Stain	Oil	Static	
Spirit Millenium	Basic	✓	✓	✓	✓	✓		250,000 Cycles #8 Duct
Neochrome II	Basic	✓	✓	✓	✓	✓		250,000 Cycles #8 Duct
Naugasoft*	Upgrade	✓	✓	✓	✓		✓	One Million Cycles #8 Duct
New Duration	Upgrade	✓	✓	✓	✓	✓	✓	400,000 Cycles #8 Duct
Royal Hampton	Basic	✓	✓		✓	✓	✓	100,000 Cycles #8 Duct
Ultrasoftouch	Upgrade		✓		✓	✓	✓	100,000 Cycles Wire Screen

FOAM QUALITY!

New High Quality Foam is a must – not an option. There is a lot of confusion in the marketplace regarding the recovering of Dental stools or chairs.

In an effort to cut costs, many upholsterers strip off your old vinyl and recover the existing foam. Good quality foam usually lasts from 10 to 15 years. So if your chair is of that age, the foam has done its time and may be breaking down; not to mention what it may have absorbed. You would be repelled by the condition of some foam stripped from old kits.

Old foam reused, or a poor grade of foam, will provide less comfort and will shorten the life and fresh appearance of your upholstery package.

Aerostatic guarantees that we use High Quality New Foam in all our upholstery kit packages, and take your old foam to the dump! End of story!